ROCK BARIATRIC CHAIR



KEY FEATURES

- → Cross-braced skis provide additional strength and manoeuvrability
- ightarrow One-piece seat with wipe-clean laminate surface
- ightarrow Strong and durable solid Birch frame
- ightarrow Seat static weight tested to 1100lbs



MATERIAL/FABRICATION

Same features as Rock dining chair just a wider seat platform making this chair ideal for bariatric use. Lacquered solid Birch frame with plywood seat in Egger Platinum White ST15 laminate finish. Key to the range's strength (tested for severe contract use; BS EN 16139:2013 level 2) are the interconnections which join each leg together, creating a rigid solid Birch frame which maximizes durability and minimizes damage.

TEST STANDARDS

- Strength and stability tested for severe contract use (BS EN 16139:2013 level 2)
- → Exceeded the 200,000 cycle pass mark by a further 300,000 cycles
- → EN1728:2012 6.4 seat and back tested to 1100lbs

OPTIONS

- → Rock Bariatric chair not available with arms
- ightarrow Matching coffee and dining tables available
- \rightarrow Matching standard size Rock dining chairs available

1RK1-BIR-WHITE 32H | 20.4W | 20.4D (in) Seat height: 18.1in Weight: 16.5lbs



COLOR/FINISH

Seat platform is in Egger Platinum White ST15 laminate finish. Frame comes in Birch as standard

ASSEMBLY

Arrives assembled and ready to use

INSTALLATION/FITTING

No fitting required

MAINTENANCE

Easy to keep clean with a soft damp cloth and use of a non-abrasive cream cleaner is acceptable. Do not use solvents, bleaches, abrasives or synthetic detergents. Avoid getting product wet. If this occurs dry thoroughly. Table tops are heat resistant but not heat proof so avoid putting items straight from the oven directly in contact with the top

WARRANTY

5 year guarantee

Pineapple

800.496.9324 inquiries@pineapplecontracts.com us.pineapplecontracts.com

E W

Т