

Pineapple

LIBERTY

Designed to enhance hygiene,
accessibility & inspection



LIBERTY



REMOVABLE SEATS
Removable seats allowing for safe inspection and thorough cleaning, securely locked in place during use

STEP 1
Insert key into back of seat and turn counter-clockwise to unlock

Our Liberty range is designed with built-in features that enhance hygiene, accessibility and inspection.



High floor clearance allows the use of hoists to manoeuvre less mobile patients



STEP 2
Push and remove back



STEP 3
Lift and remove seat



STEP 4
Clean and inspect before replacing back and seat, locking securely with key



LIBERTY




Liberty chair
 1LY1-300
 40.9"H | 33.5"W | 35.4"D
 Seat height: 18.9"
 Weight: 162 lbs.



Liberty 3-seater
 1LY3-300
 40.9"H | 81.9"W | 35.4"D
 Seat height: 18.9"
 Weight: 316.4 lbs.



Liberty 2-seater
 1LY2-300
 40.9"H | 60.6"W | 35.4"D
 Seat height: 18.9"
 Weight: 242.5 lbs.

KEY FEATURES	OPTIONS
<ul style="list-style-type: none"> → Plywood reinforcement for added strength → Removable seat allows thorough inspection and cleaning → Seat and back are securely locked in place during use 	<ul style="list-style-type: none"> → Wide choice of water-resistant upholstery 

Removable seat allows thorough inspection and cleaning



See for yourself...

To arrange a demonstration or book a visit to our showroom, please call 800.496.9324

Pineapple Contracts Inc., 2149 Avon Industrial Dr., Rochester Hills, MI 48309



800.496.9324



inquiries@pineapplecontracts.com



us.pineapplecontracts.com

Pineapple